

Thermo Fisher Scientific

En Vallaire, 1024 Ecublens, Switzerland

Tel : (41) 21 694 71 11

Field Service : (41) 21 694 75 00

Fax : (41) 21 694 71 12

E-Mail : customer@thermofisher.com

QUALITY ASSURANCE & CALIBRATION REPORT

QC/QA/CA Report

Customer : HONDA OF SOUTH CAROLINA / US

Instrument Nr. : ARL 3460 - 6788

Computer Nr. : 19011 - 18253

Place & Date :

Ecublens, the 05 March 2015

Quality assurance performed by :

Edgar Vicente

Avertissement à nos clients

Ce rapport d'essai est un document interne Thermo Fisher Scientific destiné à notre ingénieur de service. Il doit lui faciliter le travail lors de l'installation et la mise en service de notre appareil dans vos laboratoires.

Nous vous prions de bien vouloir le lui remettre à son arrivée.

Bemerkung für unsere Kunden

Dieser Versuchsrapport ist ein internes Thermo Fisher Scientific Dokument, unserm Service-Ingenieur zugehört.

Dieses soll ihm die Arbeit anlässlich der Installation sowie der Inbetriebnahme unseres Gerätes in Ihrem Labor erleichtern.

Wir bitten Sie, ihm dieses Dokument bei seiner Ankunft zu überreichen.

Advice to our customers

This testing report is an internal Thermo Fisher Scientific document intended to our service engineer. In order to make his work easier when installing and working up our instrument in your laboratories.

Please hand it over to him on his arrival.

Informations stored in OXSAS

The following instrument informations are contained in the OXSAS Database

Calibrated method(s) for concentration analysis
Methods for intensity measurements; for service engineer only
Standardisation update batch(es)
All measured values during QA/calibration process

Report of measured concentrations with standard samples

C:\Thermo\OXSAS_DATA\QA_CA_Report\C3466788.pdf

Software options list is available on the OXSAS distribution disk

Analytical programme

Instrument : ARL 3460 - 6788

Customer : HONDA OF SOUTH CAROLINA

Grating : 1080

Base(s) : Al Lines : 13

Chn	Elem	Line	Ord	Position	Slit	TRS	SDat	Al				
								Min	Max			
5	Al4	305.47	2	610.94	75			Int	Std			
1	Cr4	425.43	1	425.43	25			0.00002	1.0			
2	Cu9	510.55	1	510.55	50			0.02	15.0			
11	Fe7	371.99	2	743.98	25			0.0001	3.0			
3	Mg1	279.08	2	558.16	50			0.00005	3.0			
13	Mn5	403.45	2	806.90	25			0.00002	3.0			
4	Na1	588.99	1	588.99	37			0.000007	0.1			
9	Ni6	341.47	2	682.94	25			0.000025	5.0			
10	Pb4	368.34	2	736.68	25			0.00003	2.0			
12	Si4	390.55	2	781.10	50			0.02	27.0			
6	Sn5	317.51	2	635.02	25			0.00016	20.0			
8	Ti4	337.28	2	674.56	25			0.000006	1.0			
7	Zn4	330.26	2	660.52	37			0.00008	10.0			

Calibrations

Instrument : ARL 3460 - 6788

Customer : HONDA OF SOUTH CAROLINA

Grating : 1080

Base(s) : Al

Base : Al

x **CAL-34-AI-04** Al/Si/Cu 4xxx series

x Delivers with the certificats.

BEC & DL

Instrument : ARL 3460 - 6788

Customer : HONDA OF SOUTH CAROLINA

Grating : 1080

Base(s) : Al

Line	Wavelength		BEC_ppm	LOD_ppm	Type	Std	Window	Source Cond	Repro for LOD		RSD
	nm	Chn							Average	Sigma	
Cr0	425.43x1	1	6.37	0.108	R	Al40	NOTRS	ALARC	0.010	0.00006	0.56
Cr4	425.43x1	1	69.70	0.956	R	Al4	NOTRS	ALSPARK	0.023	0.00011	0.46
Cu9	510.55x1	2	1694.98	16.975	R	Al4	NOTRS	ALSPARK	0.023	0.00008	0.33
Fe0	371.99x2	11	13.22	0.374	R	Al40	NOTRS	ALARC	0.009	0.00009	0.94
Fe7	371.99x2	11	168.09	3.369	R	Al4	NOTRS	ALSPARK	0.024	0.00016	0.67
Mg1	279.08x2	3	14.75	0.216	R	Al4	NOTRS	ALSPARK	0.008	0.00004	0.49
Mn0	403.45x2	13	8.18	0.177	R	Al40	NOTRS	ALARC	0.011	0.00008	0.72
Mn5	403.45x2	13	77.98	2.009	R	Al4	NOTRS	ALSPARK	0.022	0.00019	0.86
Na1	588.99x1	4	4.82	0.083	R	Al4	NOTRS	ALSPARK	0.021	0.00012	0.58
Ni0	341.47x2	9	6.28	0.225	R	Al40	NOTRS	ALARC	0.007	0.00008	1.20
Ni6	341.47x2	9	55.11	1.312	R	Al4	NOTRS	ALSPARK	0.019	0.00015	0.79
Pb0	368.34x2	10	19.73	0.466	R	Al40	NOTRS	ALARC	0.007	0.00006	0.79
Pb4	368.34x2	10	223.14	2.945	R	Al4	NOTRS	ALSPARK	0.020	0.00009	0.44
Si4	390.55x2	12	1375.54	18.742	R	Al4	NOTRS	ALSPARK	0.022	0.00010	0.45
Sn0	317.51x2	6	25.90	0.699	R	Al40	NOTRS	ALARC	0.010	0.00009	0.90
Sn5	317.51x2	6	184.55	3.059	R	Al4	NOTRS	ALSPARK	0.022	0.00012	0.55
Ti0	337.28x2	8	1.97	0.035	R	Al40	NOTRS	ALARC	0.008	0.00005	0.59
Ti4	337.28x2	8	9.31	0.154	R	Al4	NOTRS	ALSPARK	0.022	0.00012	0.55
Zn0	330.26x2	7	52.29	1.082	R	Al40	NOTRS	ALARC	0.008	0.00005	0.69
Zn4	330.26x2	7	299.19	4.761	R	Al4	NOTRS	ALSPARK	0.022	0.00012	0.53

Instrument configuration

Instrument : ARL 3460 - 6788
Base(s) : Al **Lines :** 13

Customer : HONDA OF SOUTH CAROLINA

Grating : 1080

PROFIL POSITION : **Average =** 100.7 Div. **0 =** 84 Div. + 4 Turns

VACUUM PUMP : POMPE M2005SD OIL, ALCATEL (A042289)

Serial number : AM731684

VACUUM LOSS : 55 uHg/h

SOURCE + STAND : HIREP II + & MBS 301 I Manual Version V6

ARL Hirep Number S701999

Condition select board	Ponts	Setting
	W1	0
	W2	0
	W3	0
	W4	0

ARGON QUALITY : AR 48

ARGON FLOW : **Analyse Argon Flow :** 3.5 l/min

Standby Argon Flow : 0.34 l/min

Entry Pression : ~ 2.0 Bar(Kg/cm²)

FIRMWARE : **ICS :** ICS 386/U1/OE(34/44/E97) V6.0-2 14/10/04 CCS/CND-C013

CCS : NONE

SOFTWARE : OXSAS V2.1.5a

SPC Version : BASIC 9286-RP vers. 4.51.0289

Analytical configuration

Instrument : ARL 3460 - 6788

Customer : HONDA OF SOUTH CAROLINA

Grating : 1080

Base(s) : Al Lines : 13

Base	Electrode	Petrey table	Disk table	Insulator	Electrode holder	Brush
Al	3460_2-6mm S423156	3460_W_D12 S703281	None .	3460 S418820	None .	Mo S422255
Remark :						

Optical Layout material

Instrument : ARL 3460 - 6788

Customer : HONDA OF SOUTH CAROLINA

Grating : 1080

Base(s) : Al Lines : 13

Mirrors	A	B	C	D	Total
Total	2	4	3	4	13

Feet	Bot	Top	Total
Total	7	6	13

Slits	25	37	50	75	Total
Total	7	2	3	1	13

Dec.					
Total					

Filters	-10	-13	-14	-16	Total
Total	4	2	4	1	11

PMTs	39-1	39-2	40-2	Total
Total	1	9	3	13

Masks	GMask	Total
Total	2	2

Adaptors	I	II	Total
Total	5	8	13

Screws Adaptor	Clamps Adaptor
AA26818 26	S313031 13

Screws	Covers	Washers
A028013 65	S421267 65	S421267 65

Phototubes	TOP	77	73	69	65	61	57	53	49	45	41	37	33	29	25	21	17	13	9	5	1	I
		78	74	70	66	62	58	54	50	46	42	38	34	30	26	22	18	14	10	6	2	III
	BOTTOM	75	71	67	63	59	55	51	47	43	39	35	31	27	23	19	15	11	7	3	II	
		76	72	68	64	60	56	52	48	44	40	36	32	28	24	20	16	12	8	4	IV	

Optical Layout place

Instrument : ARL 3460 - 6788

Customer : HONDA OF SOUTH CAROLINA

Grating : 1080

Base(s) : Al Lines : 13

Place	Elem	Mini	Avge	Maxi	WaveLn	Ord	Positi	Mirror	Slit	Dec	Pmt	Adap	Filter	Int	Trs	Mask	Remark
1	Top		342	348	354			B				II					
2	Top		348	354	361			C				I					
3	Bot		354	361	367			B				II					
4	Bot		361	367	373			C				I					
5	Top		367	373	379			A				I					
6	Top		373	379	385			D				II					
7	Bot		379	385	391			A				I					
8	Bot		385	391	397			D				II					
9	Top		391	397	404			B				II					
10	Top		397	404	410			C				I					
11	Bot		404	410	416			B				II					
12	Bot		410	416	422			C				I					
13	Top	Cr4	416	422	428	425.43 x 1	425.43	A	25		40-2	I		1	No		
14	Top		422	428	435			D				II					
15	Bot		428	435	441			A				I					
16	Bot		435	441	447			D				II					
17	Top		441	447	453			B				II					
18	Top		447	453	459			C				I					
19	Bot		453	459	465			B				II					
20	Bot		459	465	471			C				I					
21	Top		465	471	478			A				I					
22	Top		471	478	484			D				II					
23	Bot		478	484	490			A				I					
24	Bot		484	490	496			D				II					
25	Top		490	496	502			B				II					
26	Top		496	502	508			C				I					

Optical Layout place

Instrument : ARL 3460 - 6788

Customer : HONDA OF SOUTH CAROLINA

Grating : 1080

Base(s) : Al Lines : 13

Place	Elem	Mini	Avg	Maxi	WaveLn	Ord	Positi	Mirror	Slit	Dec	Pmt	Adap	Filter	Int	Trs	Mask	Remark
27 Bot	Cu9	502	508	514	510.55	x 1	510.55	B	50		40-2	II		2	No		
28 Bot		508	514	520				C				I					
29 Top		514	520	526				A				I					
30 Top		520	526	532				D				II					
31 Bot		526	532	539				A				I					
32 Bot		532	539	545				D				II					
33 Top		539	545	551				B				II					
34 Top		545	551	557				C				I					
35 Bot	Mg1	551	557	563	279.08	x 2	558.16	B	50		39-2	II	-13	3	No		
36 Bot		557	563	569				C				I					
37 Top		563	569	575				A				I					
38 Top		569	575	581				D				II					
39 Bot		575	581	587				A				I					
40 Bot	Na1	581	587	593	588.99	x 1	588.99	D	37		40-2	II	-16	4	No		
41 Top		587	593	599				B				II					
42 Top		593	599	605				C				I					
43 Bot		599	605	611				B				II					
44 Bot	Al4	605	611	617	305.47	x 2	610.94	C	75		39-1	I	-13	5	No		Int. Std
45 Top		611	617	623				A				I					
46 Top		617	623	629				D				II					
47 Bot		623	629	635				A				I					
48 Bot	Sn5	629	635	641	317.51	x 2	635.02	D	25		39-2	II	-14	6	No		
49 Top		635	641	647				B				II					
50 Top		641	647	653				C				I					
51 Bot		647	653	659				B				II					
52 Bot	Zn4	653	659	665	330.26	x 2	660.52	C	37		39-2	I	-14	7	No	GMask	

Optical Layout place

Instrument : ARL 3460 - 6788

Customer : HONDA OF SOUTH CAROLINA

Grating : 1080

Base(s) : Al Lines : 13

Place	Elem	Mini	Avg	Maxi	WaveLn	Ord	Positi	Mirror	Slit	Dec	Pmt	Adap	Filter	Int	Trs	Mask	Remark
53 Top		659	665	671				A				I					
54 Top	Ti4	665	671	676	337.28 x 2	674.56	D	25			39-2	II	-14	8	No		
55 Bot		671	676	682				A				I					
56 Bot	Ni6	676	682	688	341.47 x 2	682.94	D	25			39-2	II	-14	9	No	GMask	
57 Top		682	688	694				B				II					
58 Top		688	694	700				C				I					
59 Bot		694	700	706				B				II					
60 Bot		700	706	712				C				I					
61 Top		706	712	717				A				I					
62 Top		712	717	723				D				II					
63 Bot		717	723	729				A				I					
64 Bot		723	729	735				D				II					
65 Top	Pb4	729	735	741	368.34 x 2	736.68	B	25			39-2	II	-10	10	No		
66 Top	Fe7	735	741	747	371.99 x 2	743.98	C	25			39-2	I	-10	11	No		
67 Bot		741	747	752				B				II					
68 Bot		747	752	758				C				I					
69 Top		752	758	764				A				I					
70 Top		758	764	769				D				II					
71 Bot		764	769	775				A				I					
72 Bot		769	775	781				D				II					
73 Top	Si4	775	781	786	390.55 x 2	781.10	B	50			39-2	II	-10	12	No		
74 Top		781	786	792				C				I					
75 Bot		786	792	798				B				II					
76 Bot		792	798	803				C				I					
77 Top	Mn5	798	803	809	403.45 x 2	806.90	A	25			39-2	I	-10	13	No		
78 Top		803	809	815				D				II					

Attenuators & Programmable Attenuators positions

Instrument : ARL 3460 - 6788

Customer : HONDA OF SOUTH CAROLINA

Grating : 1080

Base(s) : Al Lines : 13

* = Bridge cutted

' = Bottom Programmable Card

Programmable Attenuators Bit0

Programmable Attenuators Bit1

Chn	Name	Position	Prog	Base(s)	Position	Base(s)
1	Cr4	3-10				
2	Cu9	4-14				
3	Mg1	4-14				
4	Na1	2-12				
5	Al4	4-15				
6	Sn5	2-12				
7	Zn4	2-10				
8	Ti4	2-14				
9	Ni6	2-14				
10	Pb4	2-13				
11	Fe7	2-15				
12	Si4	3-13				
13	Mn5	2-15				

Table of Contents

Header	1
Advices	2
Info in OXSAS	3
Analytical programme	4
Calibrations	5
BEC & DL	6
Base AI	6
Configuration	7
Instrument	7
Analytical	8
Optical Layout	9
Material	9
Place	10
Attenuators	13